







# fortress diagnostics

An ISO 13485 accredited company

BX60792A

96 TESTS

STORE AT 2-8°C

FOR IN-VITRO DIAGNOSTICS USE ONLY

## HIV (Ag/Ab) 4th Gen HIV-1 Group M.A.B.C.D.E.F.G.H Group O, Group N & HIV-2

### ELISA

**Intended Use:**  
Fortress kit is an enzyme-linked immunosorbent assay (ELISA) for qualitative determination of antigen or antibodies to human immunodeficiency virus (HIV) type 1 and/or type 2 in human serum or plasma. The kit is intended for screening of blood donors and as an aid for the diagnosis of clinical conditions related to infection with HIV-1 and/or HIV-2, e.g., the acquired immunodeficiency syndrome (AIDS).

The human immunodeficiency viruses type 1 and type 2 are the etiologic agents of the acquired immunodeficiency syndrome (AIDS) and related conditions. HIV has been isolated from patients with AIDS, AIDS related complex (ARC) and from healthy individuals of highest risk for AIDS infection with HIV. It is followed by an acute flu-like syndrome, followed by a latent period, and then by a second phase of HIV infection. HIV infection may not be clear in many cases. The acute phase progresses to clinical AIDS in about 50% of infected individuals within 10 years after seroconversion.

Microscopic evidence of infection with HIV may be obtained by detection of virus or detection of antibodies to HIV. HIV infection is typically followed by an asymptomatic carrier state, which progresses to clinical AIDS in about 50% of infected individuals within 10 years after seroconversion.

Immunoreactive the regions of the proteins and also enabled the production of conjoined HIV-1/HIV-2 assays. The recombinant antigen could also be produced with considerably more purity and in large amounts, and they can be bound to solid-phase surfaces with much higher control over protein ratios and concentrations.

The first and second generations HIV kits were based on indirect ELISA and used only 96 antibodies only by enzyme-linked anti-human IgG antibody.

The third generation ELISA utilized double antigen "sandwich" method, again with antigen coated on solid phase polystyrene plates, but with antibodies detection achieved with the help of another enzyme-labeled antigen. The third generation assays could increase the assay's sensitivity 100 times compared to the previous generations. In addition, the detection of IgM antibodies that are present only during the early stages of infection, much shorter than the antibody detection window period (the period of time in which there is no detectable antibody production) and compare to the detection window period (days) could detect antibodies 11 days earlier.

To reduce even further the antibody detection window period, 4th generation HIV ELISAs that could simultaneously detect HIV antigens (p24) and antibodies have been developed and are commercially available. With detection of p24, the 4th generation tests shorten the window period for HIV detection to 10-15 days after infection. HIV infection could be detected 8 days earlier.

**Principle:**  
The HIV (1/2) Ag/Ab ELISA kit is a two-step incubation, "sandwich" enzyme immunoassay, which uses polystyrene microwell strips pre-coated with antigen (p24) and antibody (anti-p24) and anti-HIV (p24) and recombinant HIV-2 (gp36) and anti-HIV (p24) antibodies. As a first step, biotinylated anti-HIV (p24) antibodies together with the patient's serum or plasma sample are added into the wells. During incubation, the specific HIV-1/2 antibodies if present in the sample will bind to the antigen. In the second step, streptavidin-conjugated anti-HIV (p24) antibodies will bind to the biotinylated antibodies. The enzyme-labeled anti-HIV (p24) antibodies will bind to the streptavidin-conjugated antibodies. The enzyme-labeled antibodies will then react with the substrate to produce a color change. The amount of color intensity or T24 captured in the wells, and to the sample respectively. Wells containing samples negative for anti-HIV 1/2 or p24 remain colourless.

Component 1: 96 Tests	
HIV 1+2 (Ag/Ab) 4th Gen	Volume
Microwell Plate 96	1 plate (12x8/63/12 well strips)
Test	121ml
Negative Control	121ml
Antibody Positive Control (HIV 1)	121ml
Antibody Positive Control (HIV 2)	121ml

Antigen Positive Control	1x1ml
HIV - Conjugate Reagents	1x12ml
Biotin-Conjugate Reagent	1x3.5ml
Stock Wash Buffer	1x50ml (Dilute 1 to 20 with distilled water before use. Once diluted stable for two weeks at 2-8°C).
Chromogen Solution A	1x6ml (Ready to use and once open, stable for one month at 2-8°C).
Chromogen Solution B	1x6ml (Ready to use and once open, stable for one month at 2-8°C).
Stop Solution	1x6ml
Plastic Sealable Bag	1 Units
Plate Cover	1 Sheet
Package Inserts	1 Copy

### Additional Materials and Instruments Required But Not Provided:

1. Freshly distilled or deionized water.
2. Disposable gloves and timer.
3. Disposable waste containers for potentially contaminated materials.
4. Disposable V-shaped trough.
5. Dispensing system and/or pipette (single or multichannel).
6. Disposable tissue or clean towel.
7. Disposable gloves and clean towel.
8. Microtiter for displaying and mixing conjugate with samples.
9. Microplate reader, single wavelength 450nm or dual wavelength 450nm and 630nm.
10. Microplate agitator/shaker system.

### Sample Collection, Transportation And Storage:

1. Serum: Collect serum from the patient and store at 2-8°C. The serum should be collected by venipuncture should be allowed to clot initially and completely - the serum/plasma must be separated from the clot as early as possible as to avoid hemolysis of the RBC. Core should be taken to ensure that the serum samples are clear and not contaminated by microorganisms. Any visible particles within the sample should be removed by centrifugation at 2000g for 10 minutes. Plasma samples collected into EDTA, sodium citrate or heparin may be tested, but highly faecal, icteric, or haemolysed samples should not be used as they could give erroneous results in the assay. Do not heat or freeze samples. This can cause sample deterioration.
2. Urine: Urine samples can be used for screening. Urine samples should be collected by voiding into a clean container. Urine samples should be stored at 2-8°C. Urine samples should be tested within 3 days of collection. For shipment, samples should be packaged and labelled in accordance with the shipping local and international regulations for transport of clinical samples and biological agents.

### Special Instructions:

1. A good washing procedure is essential to obtain correct and precise analytical data.
2. It is therefore recommended to use a good quality ELISA microplate washer, manufactured at the best level of working performance. In general, no less than 3 automatic washing cycles of 350-400µl/well are required to obtain false positive reactions and high background.

3. To avoid cross-contaminations of the plate with sample or HIV-conjugates, after incubation do not discard the content of the wells but allow the plate washer to aspirate it automatically. Anyway, we recommend calibrating the washing system on the kit itself in order to match the decided analytical performance. Assume that the microplate washer liquid dispensing channels are not blocked or contaminated and sufficient volume of Wash buffer is dispensed each time into the wells.
4. In case of manual washing, we suggest to carry out 5 cycles, dispensing 350-400µl/well and aspirating the liquid for 5 times. If poor results (high background) are observed, increase the washing cycles or soaking time per well.
5. In any case, the liquid aspirated out the strips should be treated with a sodium hypochlorite solution at a final concentration of 2.5% for 24 hours, before liquid are washed in an appropriate way.
6. The concentrated washing solution should be diluted 1 to 20 before use. For one plate, mix 50ml of the concentrate with 950ml of water for a final volume of 1000ml diluted Wash Buffer. If less than a whole plate is used, prepare the proportional volume of solution.
7. Storage and Stability: The components of the kit will remain stable through the expiration date indicated on the label and package when stored between 2-8°C, do not freeze. To assure maximum performance of this HIV (1/2) Ag/Ab ELISA kit during storage, protect the reagents from contamination with microorganism or chemicals.

### Precautions And Safety:

This kit is intended FOR IN VITRO USE ONLY. [W]

### FOR PROFESSIONAL USE ONLY

The ELISA assay is a time and temperature sensitive method to avoid incorrect result, strictly follow the test procedure steps and do not modify them.

1. Do not exchange reagents from different lot, or use reagents from other commercially available kits. The components of the kit are precisely matched as to achieve optimal performance during testing.
2. Make sure that all reagents are within the validity indicated on the kit box and are of the same lot. Never use reagents beyond the expiry date stated on reagents label or on the kit box.
3. CAUTION - CRITICAL STEP: Allow the reagents and samples to stabilize to room temperature (18-20°C) before use. Samples should be gently shaken, and return to 2-8°C immediately after use.
4. Use only sufficient volume of sample as indicated in the procedure steps. Failure to do so, may cause in low intensity of the assay.
5. Do not touch the bottom exterior of the wells. Reagents or reagents may interfere with microwell reading.
6. When reading the results, ensure that the plate bottom is dry and there are no air-bubbles inside the wells.
7. Never allow the microplate wells to dry after the washing step. Immediately proceed to the next step.