

• Indicates **cobas c** systems on which reagents can be used

Order information

Creatinine Jaffé Gen.2

700 tests	Cat. No. 04810716 190
Calibrator f.a.s. (12 x 3 mL)	Cat. No. 10759350 190
Calibrator f.a.s. (12 x 3 mL, for USA)	Cat. No. 10759350 360
Precinorm U plus (10 x 3 mL)	Cat. No. 12149435 122
Precinorm U plus (10 x 3 mL, for USA)	Cat. No. 12149435 160
Precipath U plus (10 x 3 mL)	Cat. No. 12149443 122
Precipath U plus (10 x 3 mL, for USA)	Cat. No. 12149443 160
Precinorm U (20 x 5 mL)	Cat. No. 10171743 122
Precipath U (20 x 5 mL)	Cat. No. 10171778 122
Precinorm PUC (4 x 3 mL)	Cat. No. 03121313 122
Precipath PUC (4 x 3 mL)	Cat. No. 03121291 122
Diluent NaCl 9 % (50 mL)	Cat. No. 04489357 190

System-ID 07 6928 2
Code 401

Code 401

Code 300

Code 300

Code 301

Code 301

Code 301

Code 300

Code 301

Code 240

Code 241

System-ID 07 6869 3

Roche/Hitachi **cobas c** systems

cobas c 311	cobas c 501
•	•

English

System information

CREJ2: ACN 690 (Rate blanked, compensated, serum and plasma)

CRJ2U: ACN 691 (Rate blanked, urine)

SCRE2: ACN 773 (STAT, compensated, serum and plasma, reaction time: 4)

SCR2U: ACN 774 (STAT, urine, reaction time: 4)

Intended use

In vitro test for the quantitative determination of creatinine in human serum, plasma and urine on Roche/Hitachi **cobas c** systems.

Summary^{1,2,3,4,5}

Chronic kidney disease is a worldwide problem that carries a substantial risk for cardiovascular morbidity and death. Current guidelines define chronic kidney disease as kidney damage or glomerular filtration rate (GFR) less than 60 mL/min per 1.73 m² for three months or more, regardless of cause.

The assay of creatinine in serum or plasma is the most commonly used test to assess renal function. Creatinine is a break-down product of creatine phosphate in muscle, and is usually produced at a fairly constant rate by the body (depending on muscle mass). It is freely filtered by the glomeruli and, under normal conditions, is not re-absorbed by the tubules to any appreciable extent. A small but significant amount is also actively secreted.

Since a rise in blood creatinine is observed only with marked damage of the nephrons, it is not suited to detect early stage kidney disease. A considerably more sensitive test and better estimation of glomerular filtration rate (GFR) is given by the creatinine clearance test based on creatinine's concentration in urine and serum or plasma, and urine flow rate. For this test a precisely timed urine collection (usually 24 hours) and a blood sample are needed. However, since this test is prone to error due to the inconvenient collection of timed urine, mathematical attempts to estimate GFR based only on the creatinine concentration in serum or plasma have been made. Among the various approaches suggested, two have found wide recognition: that of Cockcroft and Gault and that based on the results of the MDRD trial. While the first equation was derived from data obtained with the conventional Jaffé method, a newer version of the second is usable for IDMS-traceable creatinine methods. Both are applicable for adults. In children, the Schwartz formula is used.

In addition to the diagnosis and treatment of renal disease, the monitoring of renal dialysis, creatinine measurements are used for the calculation of the fractional excretion of other urine analytes (e. g., albumin, α -amylase). Numerous methods were described for determining creatinine. Automated assays established in the routine laboratory include the Jaffé alkaline picrate method in various modifications, as well as enzymatic tests.

Test principle^{6,7,8}

This kinetic colorimetric assay is based on the Jaffé method. In alkaline solution, creatinine forms a yellow-orange complex with picrate. The rate of dye formation is proportional to the creatinine concentration in the specimen. The assay uses "rate-blanking" to minimize interference by bilirubin. To correct for non-specific reaction caused by serum/plasma pseudo-creatinine chromogens, including proteins and ketones, the results for serum or plasma are corrected by -26 μ mol/L (-0.3 mg/dL).

Creatinine + picric acid $\xrightarrow{\text{Alkaline pH}}$ yellow-orange complex

Reagents - working solutions

R1 Potassium hydroxide: 900 mmol/L; phosphate: 135 mmol/L; pH \geq 13.5; preservative; stabilizer

R2/R3 Picric acid: 38 mmol/L; pH 6.5; non reactive buffer

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Safety data sheet available for professional user on request.

Disposal of all waste material should be in accordance with local guidelines.

This kit contains components classified as follows according to the European directive 1999/45/EC.



C – Corrosive. R1 contains potassium hydroxide.

R 1: Explosive when dry. R 4: Forms very sensitive, explosive metallic compounds. R 34: Causes burns.

S 24-25: Avoid contact with skin and eyes. S 26: In case of contact with eyes,

rinse immediately with plenty of water and seek medical advice. S 35: This

material and its container must be disposed of in a safe way.

S 36/37/39: Wear suitable protective clothing, gloves and eye/face protection.

S 45: In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

Contact phone: all countries: +49-621-7590, USA: +1-800-428-2336

Reagent handling

Ready for use.

Storage and stability

CREJ2

Shelf life at 15-25 °C:

See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer:

8 weeks

Diluent NaCl 9 %

Shelf life at 2-8 °C:

See expiration date on **cobas c** pack label.

On-board in use and refrigerated on the analyzer:

12 weeks

Specimen collection and preparation⁹

For specimen collection and preparation, only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin and K₂-EDTA plasma

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Urine.

Collect urine without using additives. If urine must be collected with a preservative for other analytes, only hydrochloric acid (14 to 47 mmol/L urine, e.g. 5 mL 10 % HCl or 5 mL 30 % HCl per liter urine) or boric acid (81 mmol/L, e.g. 5 g per liter urine) may be used.

Stability in *serum/plasma*:¹⁰
 7 days at 15-25 °C
 7 days at 2-8 °C
 3 months at (-15)-(-25) °C

Stability in *urine* (without preservative):¹⁰
 2 days at 15-25 °C
 6 days at 2-8 °C
 6 months at (-15)-(-25) °C

Stability in *urine* (with preservative):¹¹
 3 days at 15-25 °C
 8 days at 2-8 °C
 3 weeks at (-15)-(-25) °C

Centrifuge samples containing precipitates before performing the assay.

Materials provided

See "Reagents - working solutions" section for reagents.

Materials required (but not provided)

See "Order information" section.

General laboratory equipment

Assay

For optimum performance of the assay, follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions. The performance of applications not validated by Roche is not warranted and must be defined by the user.

Application for serum and plasma

cobas c 311 test definition

Assay type	Rate A		
Reaction time / Assay points	10/27-37 - 15-23 (STAT 4/12-19)		
Wavelength (sub/main)	570/505 nm		
Reaction direction	Increase		
Units	μmol/L (mg/dL, mmol/L)		
Reagent pipetting	Diluent (H ₂ O)		
R1	13 μL	77 μL	
R3	17 μL	30 μL	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	10 μL	–	–
Decreased	10 μL	20 μL	80 μL
Increased	10 μL	–	–

Enter the correction value for the non-specific protein reaction as the instrument factor $y = ax + b$ for mg/dL or for μmol/L, where $a = 1.0$ and $b = -0.3$ (mg/dL) or $a = 1.0$ and $b = -26$ (μmol/L).

cobas c 501 test definition

Assay type	Rate A		
Reaction time / Assay points	10/42-52 - 24-34 (STAT 4/17-27)		
Wavelength (sub/main)	570/505 nm		
Reaction direction	Increase		
Units	μmol/L (mg/dL, mmol/L)		
Reagent pipetting	Diluent (H ₂ O)		
R1	13 μL	77 μL	
R3	17 μL	30 μL	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	10 μL	–	–
Decreased	10 μL	20 μL	80 μL
Increased	10 μL	–	–

Enter the correction value for the non-specific protein reaction as the instrument factor $y = ax + b$ for mg/dL or for μmol/L, where $a = 1.0$ and $b = -0.3$ (mg/dL) or $a = 1.0$ and $b = -26$ (μmol/L).

Application for urine

cobas c 311 test definition

Assay type	Rate A		
Reaction time / Assay points	10/27-37 - 15-23 (STAT 4/12-19)		
Wavelength (sub/main)	570/505 nm		
Reaction direction	Increase		
Units	μmol/L (mg/dL, mmol/L)		
Reagent pipetting	Diluent (H ₂ O)		
R1	13 μL	77 μL	
R3	17 μL	30 μL	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	10 μL	6 μL	144 μL
Decreased	10 μL	2 μL	180 μL
Increased	10 μL	10 μL	115 μL

cobas c 501 test definition

Assay type	Rate A		
Reaction time / Assay points	10/42-52 - 24-34 (STAT 4/17-27)		
Wavelength (sub/main)	570/505 nm		
Reaction direction	Increase		
Units	μmol/L (mg/dL, mmol/L)		
Reagent pipetting	Diluent (H ₂ O)		
R1	13 μL	77 μL	
R3	17 μL	30 μL	
Sample volumes	Sample	Sample dilution	
		Sample	Diluent (NaCl)
Normal	10 μL	6 μL	144 μL
Decreased	10 μL	2 μL	180 μL
Increased	10 μL	10 μL	115 μL

Calibration

Calibrators	S1: H ₂ O S2: C.f.a.s.
Calibration mode	Linear
Calibration frequency	2-point calibration - after reagent lot change - and as required following quality control procedures

CREJ2

Creatinine Jaffé Gen.2

Traceability: This method has been standardized against ID/MS. For the USA, this method has been standardized against a primary reference material (SRM 914).

Quality control

For quality control, use control materials as listed in the "Order information" section.

Other suitable control material can be used in addition.

Serum/plasma

For quality control use undiluted serum control material as listed above. Other suitable control material can be used in addition.

Urine

For quality control use Precinorm PUC and Precipath PUC as listed above. Other suitable control material can be used in addition.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Follow the applicable government regulations and local guidelines for quality control.

Calculation

Roche/Hitachi **cobas c** systems automatically calculate the analyte concentration of each sample.

Conversion factors: $\mu\text{mol/L} \times 0.0113 = \text{mg/dL}$

$\mu\text{mol/L} \times 0.001 = \text{mmol/L}$

Limitations – interference¹²

Criterion: Recovery within $\pm 10\%$ of initial value at a creatinine concentration of $80 \mu\text{mol/L}$ (0.90 mg/dL) in serum/plasma and $2500 \mu\text{mol/L}$ (28.3 mg/dL) in urine.

Serum/plasma

Icterus (**CREJ2**): No significant interference up to an I index of 5 for conjugated bilirubin and 10 for unconjugated bilirubin (approximate conjugated bilirubin concentration: $86 \mu\text{mol/L}$ (5 mg/dL) and approximate unconjugated bilirubin concentration: $171 \mu\text{mol/L}$ (10 mg/dL)).

Icterus (**SCRE2**): No significant interference up to an I index of 2 for conjugated bilirubin and 3 for unconjugated bilirubin (approximate conjugated bilirubin concentration: $34 \mu\text{mol/L}$ (2 mg/dL) and approximate unconjugated bilirubin concentration: $51 \mu\text{mol/L}$ (3 mg/dL)).

Hemolysis: No significant interference up to an H index of 1000 (approximate hemoglobin concentration: $621 \mu\text{mol/L}$ (1000 mg/dL)).

Lipemia (Intralipid): No significant interference up to an L index of 800. There is poor correlation between the L index (corresponds to turbidity) and triglycerides concentration.

Drugs: No interference was found at therapeutic levels using common drug panels.^{13,14}

Exception: Cefoxitin causes artificially high creatinine results.

Exception: Cyanokit (Hydroxocobalamin) may cause interference with results.

Values $< 15 \mu\text{mol/L}$ ($< 0.17 \text{ mg/dL}$) or negative results are reported in rare cases in children < 3 years and in elderly patients. In such cases use the Creatinine plus test to assay the sample.

Do not use Creatinine Jaffé for the testing of creatinine in hemolyzed samples from neonates, infants or adults with HbF levels $\geq 60 \text{ mg/dL}$ for **CREJ2** applications ($\geq 30 \text{ mg/dL}$ for **SCRE2** applications).¹⁵ In such cases, use the Creatinine plus test ($\leq 600 \text{ mg/dL}$ HbF) to assay the sample.

Estimation of the Glomerular Filtration Rate (GFR) on the basis of the Schwartz Formula can lead to an overestimation.¹⁶

In very rare cases gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.

Urine

Icterus: No significant interference up to a conjugated bilirubin concentration of $855 \mu\text{mol/L}$ (50 mg/dL).

Hemolysis: No significant interference up to a hemoglobin concentration of $621 \mu\text{mol/L}$ (1000 mg/dL).

Glucose $< 120 \text{ mmol/L}$ ($< 2162 \text{ mg/dL}$) and urobilinogen $< 676 \mu\text{mol/L}$ ($< 40 \text{ mg/dL}$) do not interfere.

Drugs: No interference was found at therapeutic levels using common drug panels.¹⁴

Exception: Cyanokit (Hydroxocobalamin) may cause interference with results.

High homogentisic acid concentrations in urine samples lead to false results.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

ACTION REQUIRED

Special Wash Programming: The use of special wash steps is mandatory when certain test combinations are run together on Roche/Hitachi **cobas c** systems. Refer to the latest version of the Carry over evasion list found with the NaOHD/SMS/Multiclean/SCCS Method Sheet and the operator manual for further instructions.

Where required, special wash/carry over evasion programming must be implemented prior to reporting results with this test.

Limits and ranges

Measuring range

Serum/plasma

$15\text{--}2200 \mu\text{mol/L}$ ($0.17\text{--}24.9 \text{ mg/dL}$)

Determine samples having higher concentrations via the rerun function.

Dilution of samples via the rerun function is a 1:5 dilution. Results from samples diluted by the rerun function are automatically multiplied by a factor of 5.

Urine

$375\text{--}55000 \mu\text{mol/L}$ ($4.2\text{--}622 \text{ mg/dL}$)

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:3.6 dilution. Results from samples diluted by the rerun function are automatically multiplied by a factor of 3.6.

Lower limits of measurement

Lower detection limit of the test

Serum/plasma

$15 \mu\text{mol/L}$ (0.17 mg/dL)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying three standard deviations above that of the lowest standard (standard 1 + 3 SD, repeatability, $n = 21$).

Urine

$375 \mu\text{mol/L}$ (4.2 mg/dL)

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying three standard deviations above that of the lowest standard (standard 1 + 3 SD, repeatability, $n = 21$).

Expected values

Serum/plasma

Adults¹⁷

Females	$44\text{--}80 \mu\text{mol/L}$	($0.50\text{--}0.90 \text{ mg/dL}$)
Males	$62\text{--}106 \mu\text{mol/L}$	($0.70\text{--}1.20 \text{ mg/dL}$)

Children¹⁸

Neonates (premature)	$25\text{--}91 \mu\text{mol/L}$	($0.29\text{--}1.04 \text{ mg/dL}$)
Neonates (full term)	$21\text{--}75 \mu\text{mol/L}$	($0.24\text{--}0.85 \text{ mg/dL}$)
2-12 m	$15\text{--}37 \mu\text{mol/L}$	($0.17\text{--}0.42 \text{ mg/dL}$)
1-< 3 y	$21\text{--}36 \mu\text{mol/L}$	($0.24\text{--}0.41 \text{ mg/dL}$)
3-< 5 y	$27\text{--}42 \mu\text{mol/L}$	($0.31\text{--}0.47 \text{ mg/dL}$)
5-< 7 y	$28\text{--}52 \mu\text{mol/L}$	($0.32\text{--}0.59 \text{ mg/dL}$)
7-< 9 y	$35\text{--}53 \mu\text{mol/L}$	($0.40\text{--}0.60 \text{ mg/dL}$)
9-< 11 y	$34\text{--}65 \mu\text{mol/L}$	($0.39\text{--}0.73 \text{ mg/dL}$)
11-< 13 y	$46\text{--}70 \mu\text{mol/L}$	($0.53\text{--}0.79 \text{ mg/dL}$)
13-< 15 y	$50\text{--}77 \mu\text{mol/L}$	($0.57\text{--}0.87 \text{ mg/dL}$)

Urine

1st morning urine¹⁷

Females	$2470\text{--}19200 \mu\text{mol/L}$	($28\text{--}217 \text{ mg/dL}$)
Males	$3450\text{--}22900 \mu\text{mol/L}$	($39\text{--}259 \text{ mg/dL}$)

Females	7000-14000 µmol/24 h	(740-1570 mg/24 h)
Males	9000-21000 µmol/24 h	(1040-2350 mg/24 h)

Creatinine clearance^{19,20} 71-151 mL/min

Refer to reference 16 for a prospective study on creatinine clearance in children.²¹

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

Precision

Precision was determined using human samples and controls in an internal protocol. *Serum/plasma*: repeatability* (n = 21), intermediate precision** (3 aliquots per run, 1 run per day, 21 days);

Urine: repeatability* (n = 21), intermediate precision** (3 aliquots per run, 1 run per day, 10 days). The following results were obtained:

Serum/plasma (CREJ2)

Repeatability*	Mean µmol/L (mg/dL)	SD µmol/L (mg/dL)	CV %
Precinorm U	105 (1.19)	2 (0.03)	2.1
Precipath U	360 (4.07)	4 (0.05)	1.1
Human serum 1	206 (2.33)	3 (0.03)	1.2
Human serum 2	422 (4.77)	5 (0.06)	1.3

Intermediate precision**	Mean µmol/L (mg/dL)	SD µmol/L (mg/dL)	CV %
Precinorm U	101 (1.14)	4 (0.05)	3.5
Precipath U	351 (3.97)	8 (0.09)	2.2
Human serum 3	201 (2.27)	5 (0.06)	2.5
Human serum 4	411 (4.64)	9 (0.10)	2.2

Urine (CRJ2U)

Repeatability*	Mean µmol/L (mg/dL)	SD µmol/L (mg/dL)	CV %
Control Level 1	8083 (91.3)	115 (1.3)	1.4
Control Level 2	15618 (177)	213 (2)	1.4
Human urine 1	19318 (218)	234 (3)	1.2
Human urine 2	7958 (89.9)	130 (1.5)	1.6

Intermediate precision**	Mean µmol/L (mg/dL)	SD µmol/L (mg/dL)	CV %
Control Level 1	8130 (91.9)	164 (1.9)	2.0
Control Level 2	15533 (176)	251 (3)	1.6
Human urine 3	19353 (219)	385 (4)	2.0
Human urine 4	7932 (89.6)	166 (1.9)	2.1

Serum/plasma (SCRE2)

Repeatability*	Mean µmol/L (mg/dL)	SD µmol/L (mg/dL)	CV %
Precinorm U	106 (1.20)	2 (0.02)	2.2
Precipath U	346 (3.91)	5 (0.06)	1.5
Human serum 1	543 (6.14)	6 (0.07)	1.1
Human serum 2	69 (0.78)	2 (0.02)	3.1

Intermediate precision**	Mean µmol/L (mg/dL)	SD µmol/L (mg/dL)	CV %
Precinorm U	100 (1.13)	4 (0.05)	4.0
Precipath U	334 (3.77)	10 (0.11)	3.0
Human serum 3	522 (5.90)	12 (0.14)	2.4
Human serum 4	64 (0.72)	3 (0.03)	5.0

Urine (SCR2U)

Repeatability*	Mean µmol/L (mg/dL)	SD µmol/L (mg/dL)	CV %
Control Level 1	6287 (71.0)	82 (0.9)	1.2
Control Level 2	15252 (172)	182 (2)	1.2
Human urine 1	24174 (273)	212 (2)	0.9
Human urine 2	2146 (24.2)	48 (0.5)	2.2

Intermediate precision**	Mean µmol/L (mg/dL)	SD µmol/L (mg/dL)	CV %
Control Level 1	6943 (78.5)	114 (1.3)	1.6
Control Level 2	15394 (174)	229 (3)	1.5
Human urine 3	24230 (274)	354 (4)	1.5
Human urine 4	2184 (24.7)	54 (0.6)	2.5

* repeatability = within-run precision

** intermediate precision = total precision / between run precision / between day precision

Method comparison

Creatinine values for human serum, plasma and urine samples obtained on a Roche/Hitachi **cobas c** 501 analyzer (y) were compared with those determined on Roche/Hitachi 917/MODULAR P analyzers (x), using the corresponding Roche/Hitachi reagent.

Serum/plasma (CREJ2)

Sample size (n) = 273

Passing/Bablok ²²	Linear regression
y = 1.000x - 0.65 µmol/L	y = 1.002x - 0.98 µmol/L
τ = 0.973	r = 0.999

The sample concentrations were between 38 and 2178 µmol/L (0.43 and 24.6 mg/dL).

Urine (CRJ2U)

Sample size (n) = 223

Passing/Bablok ²²	Linear regression
y = 0.999x + 20.66 µmol/L	y = 0.999x + 41.55 µmol/L
τ = 0.969	r = 0.999

The sample concentrations were between 934 and 50228 µmol/L (10.6 and 568 mg/dL).

Serum/plasma (SCRE2)

Sample size (n) = 224

Passing/Bablok ²²	Linear regression
y = 1.000x - 14.36 µmol/L	y = 0.996x - 12.17 µmol/L
τ = 0.964	r = 0.999

The sample concentrations were between 66 and 1775 µmol/L (0.75 and 20.1 mg/dL).

Urine (SCR2U)

Sample size (n) = 223

Passing/Bablok ²²	Linear regression
y = 0.999x + 67.83 µmol/L	y = 0.998x + 112.72 µmol/L
τ = 0.973	r = 0.999

The sample concentrations were between 931 and 48729 µmol/L (10.5 and 551 mg/dL).

References

1. Thomas C, Thomas L. Labordiagnostik von Erkrankungen der Nieren und ableitenden Harnwege. In: Thomas L, ed. Labor und Diagnose. 6th ed. Frankfurt/Main: TH-Books 2005;520-85.
2. Lamb E, Newman DJ, Price CP. Kidney function tests. In: Burtis CA, Ashwood ER, Bruns DE. Tietz textbook of clinical chemistry and molecular diagnostics. 4th ed. St. Louis, MO: Elsevier Saunders 2006;797-835.
3. <http://www.kidney.org/>
4. <http://www.nkdep.nih.gov/>
5. Lamb EJ, Tomson CRV, Roderick PJ. Estimating kidney function in adults using formulae. Ann Clin Biochem 2005;42:321-45.

CREJ2

Creatinine Jaffé Gen.2

cobas®

6. Jaffé M. Ueber den Niederschlag, welchen Pikrinsäure in normalem Harn erzeugt und über eine neue Reaktion des Kreatinins. Z Physiol Chem 1886;10:391-400.
7. Fabiny DL, Ertinghausen G. Automated reaction-rate method for determination of serum creatinine with the CentrifChem. Clin Chem 1971;17:696-700.
8. Bartels H, Böhmer M. Micro-determination of creatinine. Clin Chim Acta 1971;32:81-85.
9. Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Pre-analytical Variables. Brochure in: Samples: From the Patient to the Laboratory. Darmstadt: GIT Verlag 1996.
10. Guder WG, da Fonseca-Wollheim F, Ehret W et al. Die Qualität diagnostischer Proben, 6. Aufl. Heidelberg: BD Diagnostics, 2009.
11. Data on file at Roche Diagnostics.
12. Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
13. Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
14. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001; 38:376-385.
15. Mazzachi BC, Phillips JW, Peake MJ. Is the Jaffé creatinine assay suitable for neonates? Clin Biochem Revs 1998;19:82.
16. Filler G, Priem F, Lepage N, Sinha P, Vollmer I, Clark H, et al. β -Trace Protein, Cystatin C, β_2 -Microglobulin, and Creatinine Compared for Detecting Impaired Glomerular Filtration Rates in Children. Clin Chem 2002;48:729-736.
17. Mazzachi BC, Peake MJ, Ehrhardt V. Reference Range and Method Comparison Studies for Enzymatic and Jaffé Creatinine Assays in Plasma and Serum and Early Morning Urine. Clin Lab 2000;46:53-55.
18. Schlebusch H, Liappis N, Kalina E, Klein G. High sensitive CRP and creatinine: reference intervals from infancy to childhood. J Lab Med 2002;26:341-346.
19. Junge W, Wilke B, Halabi A, Klein G. Determination of reference intervals for serum creatinine, creatinine excretion and creatinine clearance with an enzymatic and a modified Jaffé method. Clin Chim Acta 2004;344:137-148.
20. Zawta B, Delanghe J, Taes Y, van den Noortgate N, Taes Y, Lameire N, Engel W. Arithmetic Compensation for Pseudo-Creatinine Interferences of the Creatinine Jaffé Method and its Effect on Creatinine Clearance Results. Clin Chem Part 2, Suppl S June 2001;46(6):487.
21. Wuyts B, Bernard D, van den Noortgate N, van de Walle J et al. Reevaluation of Formulas for Predicting Creatinine Clearance in Adults and Children Using Compensated Creatinine Methods. Clin Chem 2003;49:1011-1014.
22. Passing H, Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.

COBAS, COBAS C, PRECINORM and PRECIPATH are trademarks of Roche. Other brand or product names are trademarks of their respective holders. Significant additions or changes are indicated by a change bar in the margin. © 2010, Roche Diagnostics



 Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

Distribution in USA by:
Roche Diagnostics, Indianapolis, IN
US Customer Technical Support 1-800-428-2336



FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.