

Sudan University of Science and Technology
College of Graduate Studies
Department of Physics

Measurement of Temperature Using Line-width Broadening

قياس درجة الحرارة باستخدام اتساع عرض الخط الطيفي

A thesis submitted in partial fulfillment of the requirements for the degree
of M.Sc. in Physics

:Prepared by

Hassan Mansour Ahmed Daleel

:Supervised by

Dr. Mubarak Dirar Abd-Allah

ilApr 2009

Table of Contents

Page Number	Topics	No
I	Abstract	1
II	المستخلص	2
III	Acknowledgement	3
1	Chapter One: Introduction Importance of temperature 1.1 The Aim of the Work 1-2 Literature review 1-3 Presentation of the thesis 1-4	4
3	Chapter Two: Theoretical Background Prelude 2-1 Electromagnetic waves Physical properties 2-2 Atomic Spectrum 2-3 Ultra violet Units Systems 2-4 Ultra violet spectrometer 2-5 Line-width Broadening 2-6 Collision broadening 2-7 Doppler Broadening 2-8 Zeeman Splitting of Atomic of Absorbed Electromagnetic Waves 2-9	5
17	Chapter Three: Measurement of temperature Introduction 3-1 Theoretical Relation Between temperature and Doppler or 3-2 collision Broadening The proper magnetic Flux for Electromagnetic Absorption 3-3	6
20	Chapter Four: Results , Discussion and Conclusion Introduction 4-1 The Experimental set Up and Work 4-2 Results 4-3 Discussion 4-4 Conclusion 4-5 Recommendation 4-6	7
34	References	8
25	Appendix	9

Table of Figures

Page Number	Topics	No
4	The Spectrum of a Substance	2-3-1
6	Diagram of a single-beam UV/vis spectrophotometer	2-4-1
9	Block diagram of UV spectrometer	2-5-1
11	Lorentz shape and collision broadening	2-7-1
16	Energy splitting due to the effect of the magnetic field (Zeeman effect)	2-9-1
22	Relation between temperature and spectral line-width $\Delta\lambda$ for corn oil	4-3-10
22	Relation between temperature and spectral line-width $\Delta\lambda$ for sesame oil	4-3-11
23	Corn oil Spectrum at Temperature 6	4-3-1
24	Corn oil Spectrum at Temperature 12	4-3-2
25	Corn oil Spectrum at Temperature 16	4-3-3
26	Corn oil Spectrum at Temperature 22	4-3-4
27	Corn oil Spectrum at Temperature 27	4-3-5
28	Sesame Oil Spectrum at Temperature 10	4-3-6
29	Sesame Oil Spectrum at Temperature 14	4-3-7
30	Sesame Oil Spectrum at Temperature 18	4-3-8
31	Sesame Oil Spectrum at Temperature 27	4-3-9
36	Ultra Violet Spectrometer	2-5-2