APPENDICES

Appendix (1):

1. Instruments and solutions:
1.1 Instrument:
-Rotary microtome.
-Oven.
-Coplin jars.
-Staining racks.
-Microtome blade.
-Coated slides.
-Cover glass.
-Water bath.
-Dako pen.
-Thermal cycler.
-Workstation.
-Pipettors.
-Tube rack.

1.2 Solutions:

-Xylene.

-Ethyle alcohol.
-Mayers haematoxylene.
-Distilled water.
-Citrate buffer.
-Peroxidase blocker.
-Anti CDX2 antibodies (primary antibodies).
-Anti human papilloma virus antibodies (primary antibodies).
-Dextran polymer conjugated secondary antibodies and HRP.
-3.3 diaminobenzidinetetrahydrochloridin substrate buffer.
-DPX mounting media.

Appendix (3):



Cat. No. Lot No.

Quantity

2051601570 343710 7 ml

Antibody to Epstein Barr Virus (EBNA1)

Host:

Rat

Subclass:

lgG2a

Presentation:

Purified tissue culture supernatant diluted in buffer with

< 0,1 % sodium azide

Assay system:

IHC(p,f)/ELISA

Titre:

Used in an Avidin/Biotin complex AEC/DAB, it is

sufficient for 70 slides using 100µl per slide.

Fixatives:

1)NOTOXhisto 2)Formalin

Specificity:

Reacts with latent forms of EBV antigens.

Reactivity:

Human, others not tested.

Cell. Local.

Cytoplasmic

Positive control:

Epithelial carcinoma.

Storage:

Store at 2-8°C.

Errors excepted. This data sheet is a general information. The product attribute can diversify with changing Lot No as well as variations in tissue selection, tissue processing, antigen retrieval and detection systems. We do not take responsibility for any possible damage including personal injury, time or effort on economic loss caused by this product. Our warranty is limited to the price paid for the product. The product may only be used by authorized and skilled personnel. Non-application as prescribed in this data sheet leads to loss of all liability.

FOR IN VITRO USE, NOT FOR DIAGNOSTIC USE

Manufacturer

BIOCYC GmbH & Co. KG, Im Biotechnologiepark TGZ I, 14943 Luckenwalde; certified by

ISO 13485:2012 and AC: 2012 , ISO 9001:2008; Registry number: 60018916, 60018917

Distributed by:

Schichauweg 16, 12307 Berlin, Germany Tel.: ++49 (030) 765 925-0, Fax.: ++49 (030) 765 925-55 E-Mail: <u>info@quartett.com</u> Internet: www.quartett.com

Appendix (4):

Appendix 3

quartett

LOT NO. Quantity 2030415870 333313 6 ml

Antibody to CDX2

Host:

Rabbit

Immunogen:

Synthetic peptide of human CDX2 protein

Subclass:

Presentation:

Diluted antibody, in buffer with < 0.1 % sodium azide.

Assay system:

IHC(p,f)

Titre:

Used in an avidin/biotin complex AEC/DAB system vial is

sufficient for 60 slides using 100 µl per slide.

Fixatives:

1) NOTOXhisto 2)Formalin

Treatments:

Staining of formalin/paraffin tissues is enhanced by boiling the tissue sections in ProTaqs Antigen Enhancer I, Cat# 401602092 for 20 min followed by cooling at RT for 20 min.

Specificity:

CDX2 is a caudal-related homeobox transcription factort that is expressed in the nuclei of intestinal epithelium and has a function in proliferation and differentiation of intestinal cells. This antibody detects primary and metastatic colorectal carcinomas, intestinal metaplasia of the stomach

and intestinal type gastric cancer.

Reactivity:

Human, others not tested.

Storage:

Store at 2-8 °C. Do not use after expiration date indicated

on the vial.

Errors excepted. This data sheet is a general information. The product attribute can diversify with changing Lot No as well as variations in tissue selection, tissue processing, antigen retrieval and detection systems. We do not take responsibility for any possible damage including personal injury, time or effort on economic loss caused by this product. Our warranty is limited to the price paid for the product. The product may only be used by authorized and skilled personnel. Non-application as prescribed in this data sheet leads to loss of all liability.

FOR IN VITRO USE, NOT FOR DIAGNOSTIC USE

Manufacturer:
BIOCYC GmbH & Co. KG, Im Biotechnologiepark TGZ I, 14943 Luckenwalde; certified by ISO 13485:2003 and AC: 2007, ISO 9001:2008; Registry number: 60018916, 60018917

Distributed by: Jartett

Schichauweg 16, 12307 Berlin, Germany ++49 (030) 765 925-0, Fax.: ++49 (030) 765 925-55 E-Mail: info@quartett.com Internet: www.quartett.com