

## Appendix

<b>Table 2</b> <b>Section width combination for different multiple -row Detector CT configuration</b>		
Seimens	Toshiba	GE&Philips
4 x 1.25	2 x 0.5	4 x 0.5
4 x 2.5	4 x 1	4 x 1
4 x 3.75	4 x 2.5	4 x 2
4 x 5	4 x 5	4 x 3
2 x 10	2 x 8	4 x 4
2X10		4 X 5
		4 X 6
		4 X 7
		4 X 8
<p><b>Note:</b> values are number of section time section width(millimetres)</p>		