

Appendix

Table 2 Section width combination for different multiple -row Detector CT configuration		
GE&Philips	Toshiba	Seimens
4 x 1.25	2 x 0.5	4 x 0.5
4 x 2.5	4 x 1	4 x 1
4 x 3.75	4 x 2.5	4 x 2
4 x 5	4 x 5	4 x 3
2 x 10	2 x 8	4 x 4
2X10		4 X 5
		4 X 6
		4 X 7
		4 X 8
Note: values are number of section time section width(millimetres)		